
sentida hospitel

A valuable contribution to reducing the risk of falling, the Hospitel sentida low bed with the very low bedding level of 25 cm. The bedding level can be adjusted to the size of the person, even for very small residents. With a maximum working height of 80 cm, access to patients is possible even for tall nurses without the strain of bending over.







Mod.-No. 10-0462

estetica hospitel

Large double castors ensure optimum movability. The sentida Hospitel and also the estetica Hospitel are both equipped with divided side guards that support the mobilization of the resident, and thus the risk of a fall is actively prevented.




Functional elegance in nursing care

wissner-bosserhoff GmbH

Hauptstraße 4-6 D-58739 Wickede (Ruhr) Phone: +49(0)2377/784-0 Fax: +49(0)2377/784-163 E-Mail: info@wi-bo.de www.wi-bo.de



Functional elegance in nursing care





hospitel


Nursing care and hospital beds



Your advantages

at a glance:

- Can be moved and braked at any bedding height (central braking)
- Cozy design for fast recovery (healing environment)
- Mobilization-friendly side rails: without large central gap, without posts and rails
- Easy to clean thanks to plastic surfaces
- No "cross infection" through fold-down head section solution + protector-free side guard
- Automatic functions: Cardiac chair and electric legs up position
- One-handed CPR, can be reached from both sides of bed
- One-handed bed extension
- BiNetic frame for bedsore prevention and ergonomic comfort seating position:







What is hospitel?

hospital + hotel = hospitel

The hospitel concept combines comfort, convenience and multifunctionality. Its field of application is not so much acute care, but primarily the post-acute care of patients and residents in the areas of:

- Long-term care
- Geriatrics
- Orthopedics
- Rehabilitation
- Plastic Surgery

The hospitel concept supports the so-called "healing environment" and thus the rapid recovery of the patient.



High functionality

Compared to the standard range, the hospitel series of post-acute care offers the additional features:

- Fold down head-/foot sections
- Battery pack for off-grid supply
- Emergency lowering backrest (CPR) with a central head-end release
- Single-handed bed extension to 2.10 m or 2.20 m
- Horizontal and vertical round corner buffers
- Side panels with impact protection rail

