



Blues Comfort

Sofa Bed



Blues Comfort

FUNCTION

- Blues Comfort armchairs can also be used as an extra bed in patients' rooms in hospitals, hospices and so on.
- The simple unfolding mechanism made of the metal construction ensures stability when lying on the armchair.
- Wooden armrest.
- Manufactured from welded leatherette, to prevent the water to get inside the foam.



BLACK



MAHOGANY

TECHNICAL DATA

Width | 86 cm

Height | 78 cm

Depth | 98 cm

Lenght of the bed | 192 cm

Mattress platform height | 44 cm

4 Wheels for better manipulation,
2 front wheels are 7,5 cm with breaks



ALDER

WOODEN DECORS




NATURE



ALDER



CHERRY



COGNAC



TOBACCO



LARGO



WALNUT



MAHOGANY



BROWN



BLACK

COLOUR CHART



Viny 1000



Viny 6275




Viny 7095



Viny 8067



Viny 8333




Viny 6056



Viny 2120



Viny 2216




Viny 5109



Viny 5203



Viny 4057



Viny 4697



Viny 1204



Viny 1999



Viny 2450

BASIC INFORMATION

Composition | PVC 86% – PES/Cotton 14%

Weight | 610g ± 30g/m²

Width | 140 ± 3 cm

Roll length | ± 30 m

Thickness | 1,0mm ± 0,1 mm



FIRE PROOF NORM

EN 1021-1

EN 1021-2



ABRASION RESISTANCE

100,000 cycles



COLOUR FASTNESS TO ARTIFICIAL LIGHT

>_6

Blues Comfort

Sofa Bed



LINET spol. s r.o.

Želevčice 5, 274 01 Slaný, Czech Republic

tel.: +420 312 576 111, fax: +420 312 522 668, e-mail: info@linet.com, www.linet.com